

Alegria 2 *specifications*

Alegria 2 – Elevates Autoimmune and Infectious Disease Testing to the Next Level



Gain your Patient Results with high Flexibility,
Reliability and Efficiency

The new, fully automated Alegria 2 instrument provides a cost-effective sample to result solution with our unique and comprehensive Automimmune and Infectious disease panel utilising the Alegria Monotest technology.

This easy to use, customizable technique enables you to directly test any assay from the large portfolio at any time for optimized patient care.

sebia
The new language of life

System overview

Description	Fully automated immunoassay analyzer for Alegria Monotest assay processing
Parameters	Over 90 parameters available
Operation modes	Continuous or batch mode STAT samples, random access anytime
Throughput	240 tests in < 8 hours
Time to 1 st result Time to 60 results	<70 minutes <100 minutes
Sample capacity	Up to 112 primary sample tubes
Benchtop System	Integrated: - 23.6 inch TFT Full HD Touchscreen - Integrated computer - Reagents, consumables, waste container

Reagent – Consumable – Waste Management

Wash buffer container	2 liters
System fluid container	2.5 liters
Distilled water container	1 liter
Pipetting tips on board	288 tips in 3 racks
Waste capacity	Liquid waste: 10 liters Min. 180 strips Min. 250 tips
Predilutions	<i>Planned for future release</i>

Connectivity

LIS/LIMS	Bidirectional LIS LIMS (LIS2-A2 standard)
Interface	2 ethernet ports, 2 USB ports, 1 RS-232 interface
Operating system	Microsoft Windows 10 IOT
Printer	All printers supported by Windows 10 can be connected
Memory	Up to 150.000 results
System access	Individually configurable access rights and user management
Service	Integrated service software

Sample Management

Sample type	Serum - Plasma - Stool - CSF
Sample capacity	Up to 112 primary sample tubes (80 with dilution options)
Minimum sample volume	~100 µL
Sample tubes accepted	Outer diameter 9.2-16 mm Inner diameter > 8 mm Height 44-120 mm Nested microtubes
Barcodes accepted	Code 128, Code 39, Codabar/NW7, Interleave 2 of 5

Alegria Monotest Management

On-board assay portfolio	Up to 80 different assays Up to 60 strips in reaction carousel
Alegria Monotest	The trusted Alegria Monotest technology is also used on Alegria 2

Maintenance

Warm-up time	10 minutes at RT (20°C)
Close time	5 minutes: Instrument shut down and rinse
Maintenance	Daily: semi-automatic <10min Weekly: 15 minutes Monthly: 30 minutes Yearly: 1 day

Technical characteristics

Power supply	110-230Vac@60-50Hz
Power consumption	800 VA
Fuses	2 x 10.0A T (5 x 20 mm)
Dimensions	889 x 622 x 739 mm (W x D x H)
Weight	< 90 Kg / < 155 Lbs
OD wavelength	650 nm, LED
Ambient temperature	Operation: 18°C to 32°C Storage & Transportation: -20°C to 50°C
Humidity	10 to 75% without condensation
Operating altitude	Up to 3.000m
Noise level	< 65db (Decibel)

For further information, please find your Sebia local representative to contact on our Website www.sebia.com/sebia/worldwide-presence/.



ORGENTEC Diagnostika GmbH | Carl-Zeiss-Straße 49-51 | 55129 Mainz | Deutschland | Tel. : +49 6131 / 9258-0 |
Fax : +49 6131 / 9258-58 | E-mail : orgentec@orgentec.com | www.orgentec.com

©Sebia 2025 **ORGENTEC Diagnostika GmbH** - Carl-Zeiss-Straße 49-51 | 55129 Mainz | Deutschland - Ref. : QM210178 - A2FIEN - 01/2025 - 05 - Non-contractual illustration
ORGENTEC reserves the right to modify at any time and without notice the information contained in this document intended for healthcare professionals. This document contains information intended for wide distribution and may therefore contain product details or information that is not available or valid in your country. Carefully read the instructions in the reagent package inserts and instrument manuals.